

MUA 1.6 ACCESSORIES

FOR ALL SERIES 300, 400 & 500 MUA 1.6

NOTE: The MUA 1.6 system is only compatible with the MUA 1.6 accessories on this page. Older MUA accessories are incompatible with MUA 1.6.



Healing Cap

Placed on top of Multi-Unit Abutments during the gingival healing period.

MUA-0345-HCM



Placement Tool

Assists placement of all angled Multi-Unit Abutments.

MUA-0345-PTM



Open Tray Impression Coping

IPO-0345-00M



Closed Tray Impression Coping

IPC-0345-00M



Analog

Used in the dental laboratory to represent the MUA connection. New design includes a gingival mask retention groove.

ALG-0345-00M



Digital Scan Body

Snap onto the MUA for digital optical scans. Made of PEEK Medical Grade Thermoplastic.

DEN-0345-SB



Ti-base

Used by labs to make screw-retained connections to MUA16.

MUA-0345-TBM



Verification Cylinder

Used for temporary restorations and verification jigs.

MUA-0345-VCM



1.6mm Prosthetic Screw

Replacement screw for attaching Ti-base to MUA 1.6.

THS-1635-23



Accessory Torque Drivers

#6 Hexalobular Ball

REQUIRED
FOR THESE
ACCESSORIES



TDB-0060-24 TDB-0060-28