

MUA ACCESSORIES

FOR ALL SERIES 300, 400 & 500 MULTI-UNIT ABUTMENTS



Healing Cap

Placed on top of Multi-Unit Abutments during the gingival healing period.

MUA-0345-HC



Placement Tool

Assists placement of all angled Multi-Unit Abutments.

MUA-0345-PT



Impression Coping

Open tray impression post.

MUA-0345-IP



Analog

Used in the dental laboratory to represent the MUA connection.

ALG-0345-00



Digital Scan Body

Attaches to the MUA for digital optical scans.

MUA-0345-SB



Ti-base

Used by labs to make screw-retained connections to MUA.

MUA-0345-TB



Burnout Sleeve

For analog prothesis production.

MUA-0345-BN



Verification Cylinder

Used for temporary restorations and verification jigs.

MUA-0345-VC



REQUIRED FOR THESE ACCESSORIES



TDB-0060-24 TDB-0060-28

Accessory Torque Drivers

#6 Hexalobular Ball