



MULTI-UNIT ABUTMENT (MUA), STRAIGHT

Used to secure multi-unit screw-retained prosthetics. Straight MUAs are self-retaining with the threaded shank on the abutment.



PLATFORM HEIGHT	300 series	400 series	500 series
0.5 mm	MUA-0305-00	MUA-0405-00	MUA-0505-00
1 mm	MUA-0310-00	MUA-0410-00	MUA-0510-00
2 mm	MUA-0320-00	MUA-0420-00	MUA-0520-00
3 mm	MUA-0330-00	MUA-0430-00	MUA-0530-00
4 mm	MUA-0340-00	MUA-0440-00	MUA-0540-00
5 mm	MUA-0350-00	MUA-0450-00	MUA-0550-00




Use the DEH-0910-32 MUA Driver to install.



Anodization for improved aesthetics.



Series color coding.



Specifications
Material: 6AL4V E.L.I. Titanium
Profile Diameter: 3.75 mm
Platform Heights (mm): .5, 1, 2, 3, 4, 5
Torque: 20 Ncm (300 Series)
30 Ncm (400 & 500 Series)

MUA Accessories

See page 22 for details.

