


## POSITIONER DENTURE RETAINING ABUTMENT

Abutments for retaining overdenture restorations. Retention inserts sold separately.



PLATFORM HEIGHT	300 series	400 series	500 series
0.5 mm	ALR-0305-00	ALR-0405-00	ALR-0505-00
1 mm	ALR-0310-00	ALR-0410-00	ALR-0510-00
2 mm	ALR-0320-00	ALR-0420-00	ALR-0520-00
3 mm	ALR-0330-00	ALR-0430-00	ALR-0530-00
4 mm	ALR-0340-00	ALR-0440-00	ALR-0540-00
5 mm	ALR-0350-00	ALR-0450-00	ALR-0550-00



Gold anodization for superior aesthetics.



Use the Positioner Driver to attach to the implant.

**DLR-0345-00**

### Specifications

Material: 6AL4V E.L.I. Titanium  
 Profile Diameter: 3.75 mm  
 Platform heights (mm): .5, 1, 2, 3, 4, 5  
 Maximum torque: 20Ncm 300 Series  
 30Ncm 400 & 500 Series

