


RESTORATIVE KIT

KIT-REIN-00




20 Ncm Torque Wrench for Series 300 Abutments
TWB-3000-20



Positioner Denture Retaining Abutment Driver DLR-0345-00



MUA Driver DEH-0910-32



30 Ncm Torque Wrench for Series 400 & 500 Abutments
TWB-4500-30



Knurled Tapered .050 Driver
HDT-5050-02



Ball Abutment Driver
DEH-0000-00

Short Tapered .050 Drivers



HDT-0050-S0



HDT-0050-S1

Hexalobular Drivers for MUA attachments



TDB-0060-24



TDB-0060-28

Long Tapered .050 Drivers



HDT-0050-01



HDT-0050-00

Abutment Extractors



EXT-0080-00 for Series 300 abutments



EXT-0172-00 for Series 400 & 500 abutments

