


## TEMPORARY ABUTMENT

Temporary abutments are designed to be temporarily seated with temporary attachment or crown while the final restoration is being made.



PLATFORM HEIGHT	300 series	400 series	500 series
0.5 mm	ATM-0305-00	ATM-0405-00	ATM-0505-00
1 mm	ATM-0310-00	ATM-0410-00	ATM-0510-00
2 mm	ATM-0320-00	ATM-0420-00	ATM-0520-00
3 mm	ATM-0330-00	ATM-0430-00	ATM-0530-00
4 mm	ATM-0340-00	ATM-0440-00	ATM-0540-00
5 mm	ATM-0350-00	ATM-0450-00	ATM-0550-00



### Specifications

Material: 6AL4V E.L.I. Titanium  
 Profile Diameters (mm): 4, 5, 5.5  
 Platform heights (mm): .5, 1, 2, 3, 4, 5