



IMPLANT ONE **TI-BASE ABUTMENT**

Implant One Ti-Base abutments are compatible with popular CAD/CAM systems and restoration materials.

- Made from 6AL4V E.L.I. Titanium.
- Titanium retaining screw included.
- Compatible with exocad, 3Shape, and Sirona.



USAGE



Specifications
Material: 6AL4V E.L.I. Titanium
Profile Diameter: 4.5 mm
Platform height: .5 mm – 2.5 mm

Gold anodization helps improve the aesthetics on ceramic restorations

Color-coded anodization for easy platform identification.



Sizes Available and Part Numbers

PLATFORM HEIGHT	300 series	400 series	500 series
0.5 mm	ATB-0305-00	ATB-0405-00	APT-0505-00
1.5 mm	ATB-0310-00	ATB-0410-00	APT-0510-00
2.5 mm	ATB-0320-00	ATB-0420-00	APT-0520-00