



WIDE POST ABUTMENT

Wide Post Abutments are designed for being manually customized by dental laboratories for patient specific abutments.



PLATFORM HEIGHT	300 series	400 series	500 series
0.5 mm	AWP-0305-00	AWP-0405-00	AWP-0505-00
1 mm	AWP-0310-00	AWP-0410-00	AWP-0510-00
2 mm	AWP-0320-00	AWP-0420-00	AWP-0520-00
3 mm	AWP-0330-00	AWP-0430-00	AWP-0530-00
4 mm	AWP-0340-00	AWP-0440-00	AWP-0540-00
5 mm	AWP-0350-00	AWP-0450-00	AWP-0550-00
PROFILE DIAMETER	4.5 mm	5 mm	5.5 mm

Specifications

Material: 6AL4V E.L.I. Titanium
 Profile Diameters: 4.5 mm, 5 mm, 5.5 mm
 Platform heights (mm): .5, 1, 2, 3, 4, 5