



IMPLANT

— O N E —

PRECISION DENTAL IMPLANTS

SERIES

300

400

500

www.implantlogistics.com

**Provide your patients the top dental
implant system on the market today**

IMPLANT ONE MAIN FEATURES

Morse Taper Connection

Super-tight taper design dramatically reduces micro movement and micro gaps at the implant-abutment connection, creating a bacteria-free environment.

Sub-Crestal Placement

The bacteria-free connection of Implant One implants allows clinicians to place them sub-crestally, which reduces the stress applied to the cortical plate, and reduces tissue and bone die-back around the implant-abutment connection.

External Threads Designed to Reduce Torque

The Implant One system provides clinicians with two thread variants to attain optimum primary stability in any clinical case. Self-tapping threads reduce the amount of torque required for placement, which means healthier bone, and happier patients.

Complete Suite of Tools

The Implant One system includes fully-integrated surgical and restorative tools and tool kits designed to give doctors everything they need to successfully place our implants, abutments, and accessories with minimal bone and tissue loss.



IMPLANT

— O N E —

PRECISION DENTAL IMPLANTS



WIDE THREAD

The Wide Thread is ideal for immediate implant placement in extraction sites or in cases with poor bone density. The larger surface area of the threads provide maximum primary stability.

STANDARD THREAD

Suitable for all bone density types, the Standard Thread implants feature aggressive, self-tapping threads for increased initial stability and "Microthreads" to reduce pressure on the thin cortical plates, lessening the risk of pressure necrosis.

SOLUTIONS FOR EVERY CHALLENGE



ABUTMENTS

Implant One offers a complete line of abutments for individual restorations, and full-arch fixed or removable options.

Compatible with:
3shape[®] exocad

TOOLS AND KITS

Every tool you need to install and restore Implant One implants and abutments.



Implant One Surgical Kit, one of four Implant One Tool Kits




3-FLUTED DRILLS

Engineered for efficiency and bone health

ABUTMENT EXTRACTOR TOOL


Abutments are easily retrievable using the Abutment Extractor Tool included in the Restorative Kit.



SUPERIOR SCAN BODY

Top half designed for super-accurate optical scans.

Bottom half designed for radiographic confirmation of proper seating.



A PRECISION DENTAL IMPLANT SYSTEM



IMPLANT
— ONE —
PRECISION DENTAL IMPLANTS

The Implant One system
drastically reduces bone
resorption post-placement
and virtually does away
with bacteria-based issues.

ABOUT



IMPLANT
LOGISTICS
INNOVATIVE DENTAL TECHNOLOGY

**MAKERS OF THE IMPLANT ONE
IMPLANT SYSTEM**

Implant Logistics was founded by general dentist and implantologist, Dr. Leo Malin, in response to the need he saw for affordable US-based implant systems that are simple to use, preserve bone, and provide the full range of restorative options.



Dr. Leo Malin (right) and co-owner Tom Arendt

608-498-4855

sales@implantlogistics.com