




IMPLANT ONE STANDARD STRAIGHT ABUTMENT

Straight abutments are ideal for placement in narrow spacing areas.

- Made from 6AL4V E.L.I. Titanium.
- The top of the abutment is thinner than on the standard abutment.
- Straight abutments can be prepped and fitted with a custom crown.




Color-coded anodization for easy platform identification



Specifications

Material: 6AL4V E.L.I. Titanium
 Profile Diameter: 3.1 mm – 4.5 mm
 Height from top of implant to top of abutment: 8 mm



Part Numbers:

300 series	400 series	500 series
AST-0300-00	AST-0400-00	AST-0500-00
20 Ncm	25 Ncm	25 Ncm

Max. Torque:

