


## IMPLANT ONE **TEMPORARY ABUTMENT**

Designed to be temporarily seated with a temporary attachment or crown while the final restoration is being made.


- Made from 6AL4V E.L.I. Titanium.
- Available in 6 different cuff heights ranging from .5 mm to 5 mm.
- Titanium retaining screw included.
- Temporary crown must be configured to prevent loading.



Barrel can be modified chair-side. Sterilization must be done after any modification.

**Temporary abutments must never be used for permanent restoration as the connection is not sufficiently tight to prevent bacterial intrusion.**

Special connection profile prevents abutment from seating permanently in the implant.



Specifications
Material: 6AL4V E.L.I. Titanium
Profile Diameter: 4 mm – 5.5 mm
Cuff height: .5 mm – 5 mm

### Sizes Available and Part Numbers

CUFF HEIGHT	300 series	400 series	500 series
0.5 mm	ATM-0305-00	ATM-0405-00	ATM-0505-00
1 mm	ATM-0310-00	ATM-0410-00	ATM-0510-00
2 mm	ATM-0320-00	ATM-0420-00	ATM-0520-00
3 mm	ATM-0330-00	ATM-0430-00	ATM-0530-00
4 mm	ATM-0340-00	ATM-0440-00	ATM-0540-00
5 mm	ATM-0350-00	ATM-0450-00	ATM-0550-00
PROFILE DIAMETER	4 mm	5 mm	5.5 mm